

<i>Product code</i>	<i>Glass plate size (cm)</i>	<i>Number of gels that can be run per unit</i>	<i>Gel thickness (mm)</i>	<i>Number of samples/gel</i>	<i>Heat exchanger included</i>	<i>Compatible caster</i>
SE660	18 x 24 18 x 16 18 x 8	1 to 4	0.75, 1.0, 1.5	Up to 28	Yes	SE6015
SE600X	18 x 16	1 to 4	0.75, 1.0, 1.5	Up to 28	Yes	SE6015, SE615, SE675
SE600	18 x 16	1 to 4	0.75, 1.0, 1.5	Up to 28	Yes	SE6015, SE615, SE675
SE640	18 x 8	1 to 4	0.75, 1.0, 1.5	Up to 28	No	SE6015
SE410	18 x 24	1 to 2	0.75, 1.0, 1.5	Up to 28	No	SE6015
SE400	18 x 16	1 to 2	0.75, 1.0, 1.5	Up to 28	No	SE6015, SE615, SE675
SE300 miniVE	10 x 10.5 10 x 8	1 to 2	0.75, 1.0, 1.5	Up to 18	No	SE245, SE235
SE250	10 x 8	1 to 2	0.75, 1.0, 1.5	Up to 18	Yes	SE215, SE245, SE275
SE260	10 x 10.5 10 x 8	1 to 2	0.75, 1.0, 1.5	Up to 18	Yes	SE245, SE235